

STANDARD OPERATING PROCEDURE (SOP)

CLEANING OF NEAR PATIENT HEALTHCARE EQUIPMENT



Appendix 7 - Resuscitation Checklist

HOSPITAL:								
EMERGENCY DEPARTMENT								
Resuscitation Area:								
Patient trolley monitoring and work surfaces are cleaned after every patient								
Date and Time:								
Signature:								
Clean work surfaces								
Clean monitors and computer								
Clean patient chairs								
Clean procedure trolleys								
If trolley remains in the department inspect and clean * mattress								
Clean trolley frame and bed rails								
Check hand rub nozzle and replace if required								
Clean overhead lamp and check it is working								
If used clean otoscope/ auroscope								
Re-stock cupboards								
Re-make trolley with clean linen if not stored in a public area								
If used clean Bair Huggers								
If used clean Ultrasound								
If used clean Fluid Warmers								
If used clean IV Stands and Pumps								
PLEASE ENSURE ANY FAULTY EQUIPMENT IS REPORTED TO NURSE IN CHARGE AS SOON AS IT IS DISCOVERED								

* If a mattress is found to be contaminated, remove from use and inform the SCN. The SCN will inform the Lead Nurse or Service Manager.

The most up-to-date version of this policy can be viewed at the following website:

www.nhsggc.org.uk/infectioncontrol