

How to use your Genuair®

You have been prescribed a Genuair® to help manage your lung condition.

1. To open the Genuair®, lightly squeeze the arrows marked on each side of the mouthpiece cover and pull off.
2. Hold the Genuair® with the mouthpiece towards you and the coloured button facing straight up. Do not tilt.
3. Press the coloured button all the way down and then release. Do not continue to hold the coloured button down.



4. Check that the coloured window at the front of the Genuair® is green. The dose is now ready.
5. Breathe out as far as is comfortable, away from the Genuair® device.
6. Put your lips tightly around the mouthpiece. Breathe in as strongly and deeply as you can until your lungs are full.
7. You will hear a click to tell you the dose is released. Keep breathing in after you have heard the click.
8. Remove the Genuair® and hold your breath for 5 to 10 seconds and then breathe out slowly.
9. Breathe normally.
10. Make sure the coloured window has turned red. This confirms you have inhaled the dose correctly. If still green, repeat steps 4 – 9.
11. To close, replace the mouthpiece cover.

Always:

- Check the dose counter on the front of the Genuair®. This tells you approximately how many doses are left. The dose counter counts down slowly in intervals of 10, i.e. (60, 50, 40, 30, 20, 10). The last ten doses appear with a red-striped band.
- Make sure that you order a spare Genuair® before your present inhaler runs out.
- Make sure that the Genuair® is in date. The expiry date is on the Genuair® label.
- Replace after 60 days and return any partially used inhalers to your pharmacy.
- Take your inhalers with you when you attend your GP, practice nurse or the hospital.
- Keep your Genuair® dry and clean.

Cleaning:

- If necessary, wipe the mouthpiece with a dry single-use tissue. Do not wash with water.

If you need further help with your inhaler please contact your doctor, practice nurse, pharmacist or healthcare professional who gave you this leaflet.

Please use this information with the manufacturer's instruction leaflet.